

## Tufts Office of Sustainability

### ULT Freezers Best Practices Checklist

Action	Rationale
<p><b>1.1 – Keep freezers full for maximum energy efficiency</b></p> <ul style="list-style-type: none"> <li>• Consolidate freezer space</li> <li>• Share freezers with other labs</li> <li>• If there is unused space, put in reusable ice packs</li> </ul>	<ul style="list-style-type: none"> <li>• One -80°C freezer uses as much energy as a single family home (18-28 kWh/day)</li> <li>• A fully loaded freezer has minimum temperature variation</li> </ul>
<p><b>1.2 – Set ULT freezers at the highest required temperature</b></p> <ul style="list-style-type: none"> <li>• Only cool samples to -80°C if absolutely needed</li> <li>• Set ULT freezers at -70°C</li> </ul>	<ul style="list-style-type: none"> <li>• Raising the temperature by 10 degrees saves the power used by the equivalent of one standard (-20°C) freezer</li> <li>• -70°C is adequate for most bio-molecules and many microbial cultures; DNAs can be store at -20°C*</li> </ul>
<p><b>1.3 – Inventory freezer contents</b></p> <ul style="list-style-type: none"> <li>• Clearly label the samples with the: date, name of sample, researcher</li> <li>• Attach an inventory to the outside of the door that includes the location of samples</li> </ul>	<ul style="list-style-type: none"> <li>• Minimize the time the door is open while assessing samples</li> <li>• Every time a ULT freezer door is opened it takes 10 min to recover its temperature back to its set point</li> </ul>
<p><b>1.4 – Defrost regularly</b></p> <ul style="list-style-type: none"> <li>• Defrost freezers at least annually, or when the frost prevents the door from sealing properly, or when ice build-up exceeds 1 inch within the cabinet near samples.</li> </ul>	<ul style="list-style-type: none"> <li>• Clear space for sample storage</li> <li>• If frost prevents proper door seal, it increases energy consumption and affects freezer performance</li> <li>• Ice build-up may block the removal of heat</li> </ul>

<p><b>1.5 – Discard unneeded samples</b></p> <ul style="list-style-type: none"> <li>• Clean out old unwanted samples at least annually</li> </ul>	<ul style="list-style-type: none"> <li>• Normally at least 30% of the samples stored in ULFs are out-of-date or not useful for research</li> <li>• Free up space and avoid the need to buy new freezers</li> </ul>
<p><b>1.6 – Place freezers in a well-ventilated place</b></p> <ul style="list-style-type: none"> <li>• Avoid sources of heat such as sunlight or warm rooms</li> <li>• Leave 8” clear space on the top and around ULT freezers; Do not store items on top or around freezers</li> </ul>	<ul style="list-style-type: none"> <li>• Operating ULT freezers in warm environments (&gt;85°F) increases energy consumption by as much as 24%</li> <li>• Poor heat dissipation is the leading cause of component failure leading to downtime and high repair bills</li> </ul>
<p><b>1.7 – Clean dusty condenser filters</b></p>	<ul style="list-style-type: none"> <li>• Freezers with dusty condenser filters use more energy to remove heat</li> </ul>
<p><b>1.8 – Purchase efficient ULT freezers to replace old dysfunctional ones</b></p> <ul style="list-style-type: none"> <li>• In Boston campus, new units that use under 25 kWh/day may be eligible for a utility rebate. Contact Office of Sustainability for more information on the rebate.</li> </ul>	<ul style="list-style-type: none"> <li>• When new, conventional ULT freezers use about 16 to 22 kWh per day. After years of service, many consume over 30 kWh per day</li> <li>• Energy efficient units use as low as 8 kWh/day</li> </ul>
<p><b>1.9 – Install alarm device to monitor failure, if current units don’t have it</b></p> <ul style="list-style-type: none"> <li>• New freezers should have built-in alarm system</li> <li>• Install ULT freezer monitors for older units</li> </ul>	<ul style="list-style-type: none"> <li>• Monitor freezer performance and prevent failures that threaten your sample safety</li> </ul>

\*Click the link below to see a list of samples currently stored safely at -70°C across the US.

<https://www.freezerchallenge.org/score-sheet--resources.html>